EN Whitepost associated with a multifunction resin cement to simplify dental rehabilitation

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QUEJA

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EN CLINICAL/RADIOGRAPHIC EXAM

During the exam, a ceramic veneer with color alteration was observed on tooth 11 as well as cervical and incisal unevenness when compared to tooth 21. Radiographically, tooth 11 had already received a satisfactory endodontic treatment.

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TREATMENT

The treatment was initiated with the increase of the clinical crown of tooth 11 to align the margin and the gingival zenith. After tissue maturation, the pre-existing restoration was removed and a fiberglass intraradicular retainer (Whitepost DC) was installed by using the direct technique. After cementation of the post and core build up with Allcem Core, a ceramic total crown was installed.

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7 to 9. Conduct preparation with the drill available in the system (Whitepost DC #1), providing an anatomic preparation equivalent to the post format.

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Whitepost

APPROVED BY THE MOST EXIGENT PROFESSIONALS

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DOUBLE TAPER
DOBLE CONICIDAD

White Post fits better into to the conduit with no need for excessive wear of the conduit walls.

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HIGHLY TRANSLUCENT
EXCELENTE TRANSLUCIDEZ

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Photo of the drill magnified 25x.
Foto de la fresa aumentada 25x.

FGM



10 to 12. Fiberglass post adaptation into the conduct.

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
13. Application of 37% phosphoric acid (Condac 37) on the crown enamel for 30 seconds followed by thoroughly washing.

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
14. Drying with absorbing paper.

14. Secado con papel absorbente.

15. Ambar was actively applied with Cavibrush on the crown portion as well as intra-conduct with light air jets to evaporate the solvent and intra-conduct excess removal with absorbing paper.

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17. The post was treated with alcohol previously to silane (Prosil) application on the surface, waiting for 60 seconds for the solvent to evaporate.

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18. Allcem Core selection.

18. Selección del cemento AllCem Core.

19. After the adhesive preparation, insertion of the cement into the conduct with the aid of an applicator tip.

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20. Post insertion.

20. Inserción del poste.

21. Light-curing for 60 seconds, after the addition of the material for the core build up.

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22. Cemented post ready to be cut under refrigeration.

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
23. Preparation of the teeth for full crown.

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

25 and 26. Molding.

25 y 26. Moldeo.


27 and 28. With the crown already made, the adhesive preparation is made by applying 37% phosphoric acid (Condac 37) for 15 seconds and the adhesive (Ambar) in two layers and light curing for 20 seconds.

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29. In the piece interior, acid etching is performed for 20 seconds with Condac Porcelana 5% and, after that, washing and drying are carried out before silane application (Prosil). Once the adhesive preparation is done, the dual-curing resin cement (Allcem Core) is applied in the interior of the piece which is positioned on the tooth.

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30 and 31. After the excess removal, additional light-curing on the margins and finishing on the cervical terminal, one can notice a well-adapted crown and an excellent result.

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FIBERGLASS POST
Whitepost

POSTE DE FIBRA DE VIDRIO
Whitepost

5th GENERATION LIGHT-CURING ADHESIVE
Ambar

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DUAL-CURING RESIN CEMENT
Allcem Core

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HYDROFLUORIC ACID AT 5%
Condac Porcelana 5%

ÁCIDO HIDROFLUORÍDRICO AL 5%
Condac Porcelana 5%

SILANE
Prosil

SILANO
Prosil