FGM

CE 2797

Allcem cimento DUAL

BR | ES | EN | FR

PROFISSIONAL
PROFESIONAL
PROFESSIONAL
PROFESSIONNEL

DENTSCARE LTDA - Av. Edgar Nelson Meister, 474 - Distrito Industrial - 89219-501 - Joinville - SC - Brasil / Brazil. - CNPJ / Tax ID: 05.106.945/0001-06 INDÚSTRIA BRASILEIRA / BRAZILIAN INDUSTRY - Responsável Técnico / Technical Authority: Friedrich Georg Mittelstädt CRQ: 13100147 - SC - ANVISA nº 80172310032 - Marca / Brand: FGM® / Fabricado no Brasil / Made in Brazil

REP NUNO FLORES - Al. Bonifacio Lazaro Lozano, 3 – Piso 0 – C 2780-125 Oeiras / Portugal · Tel: (351) 21 4439292 · ec.rep@fgm.ind.br

Serviço ao profissional (Brasil)
0800 644 6100
Technical support (Europe)
Servicio al profesional (Americas)
+55 47 3441 6100
www.fgm.ind.br | fgm@fgm.ind.br

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Cimento resinoso dual

Somente uso profissional

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Dual-curing resin cement

Professional use only

Read all the information in this document carefully before using the product. Keep the document for later reference at least until the product is finished or until there is no more interaction between the product and its last patient.

Product description

Allcem is a dual-curing radiopaque permanent resin cement system indicated for the cementation of indirect restorations to the tooth structure. The product is composed of a base paste and a catalyst paste, stored in a double-body syringe, which allows releasing the correct proportion of each paste (1:1) and separating syringes of base paste and catalyst paste. Allcem is composed of methacrylate monomers, such as TEGDMA and BisEMA, inorganic load, photoinitiators, coinitiators, catalysts and pigments. Allcem presents adjusted viscosity allowing the formation of fine layers of cement with the appropriate covering of the parts.

Presentation forms

· Dual-body syringe (base + catalyst) with 5g and 5 self-mixing tips.
Available shades: A1, A2, A3 (Universal) and Trans.
· Simple syringe of base paste with 2.5 g.
Available shades: A1, A2, A3 (Universal) and Trans.
· Simple syringe of catalyst paste with 2.5 g.
· Packaging with 20 disposable self-mixing tips.


Basic composition

· Base paste: methacrylate monomers such as BisGMA, BisEMA and TEGDMA, camphorquinone, coinitiators, barium-aluminum-silicate glass microparticles, particles of silicon dioxide, inorganic pigments and preservatives.
· Catalyst paste: methacrylate monomers, dibenzoyl peroxide and stabilizers aluminum-silicate glass microparticles. The mixture of base paste and catalyst pastes at proportion (1:1) presents 66 - 67% load in weight.

Product indication

Allcem is a dual-curing radiopaque permanent resin cement system indicated for the cementation of indirect restorations to the tooth structure. It can be used for permanent cementation of;
· Porcelain fused-to-metal crowns and bridges;
· Metal crowns, bridges, inlays and onlays (high noble, noble, and base metals);
· Crowns and bridges with minimal tooth structure;
· Maryland bridges (resin bonded bridges);
· All ceramic/porcelain and pre-cured composite crowns, bridges, inlays and onlays;
· Endodontic posts;
Note: For cementation of laboratory-fabricated veneers, use Allcem Veneer APS.

Warnings and Precautions
During product handling
· Allcem cement contains acrylate resins. Avoid use of this product on patients with known acrylate allergies. To reduce the risk of allergic response, minimize exposure to these materials. In particular, avoid exposure to uncured resins;· Use of protective gloves and a no-touch technique is recommended;
· If skin contact occurs, wash skin with soap and water;
· Acrylates may penetrate commonly-used gloves. If cement contacts glove, remove and discard glove, wash hands immediately with soap and water and then re-glove;
· If accidental contact with eyes or prolonged contact with oral soft tissue occurs, flush with large amounts of water. If irritation persists, consult a physician.
During product application
· If there is pulp exposure, protect the area with a base/liner material;
· Do not use eugenol-based provisional cements before using Allcem;
· Allcem is not compatible with single-step self-etching adhesives, except Ambar Universal APS (FGM), which can be used according to its protocol;
· Light cure the resin cement for at least 40 seconds on each margin, with light curing unit (halogen light or LED) with power of 450mW/cm^2 or more and blue light emission spectrum of (400-500 nm). Wear protective goggles when operating a light-curing unit;
· Discard the self-mixing tip after use;
· Be sure to place back the cap after use.
Contraindication
The product should not be used in patients allergic to any of the substances contained in the formula.
Side effects
Systemic side effects are not known to date. In allergic individuals, allergic reactions may occur.
Instructions for use
Coupling the self-mixing tip to the double-body syringe:

a) After removing the cap of the double-body syringe, align the triangular-shaped protruding mark at the base of the self-mixing tip with the mark at the top of the syringe. The alignment and coupling of the two parts must occur in this position.
b) After coupling, spin (clockwise) the self-mixing tip ¼ of a turn for the correct final coupling. Make sure the coupling is perfect and ready for use.
Obs.: The self-mixing tip with applicator nozzle is disposable (use only once).
Tooth Preparation
1. Isolate the tooth and use a retraction system when the margin is subgingival to guarantee a preparation free of contamination.
2. Clean the preparations appropriately with pumice stone and prophylaxis cup, etch enamel and dentin for 15 seconds with phosphoric acid at 37% (Condac 37, FGM). Rinse with water abundantly for 10 seconds. After making sure that the whole acid solution has been removed, dry the preparation slightly with air leaving the tissues moist.
3. Apply the adhesive according to the manufacturer's instructions and light cure it. The dual-curing resin cement is not compatible with single-step self-etching adhesives, except Ambar Universal APS (FGM), which can be used according to its protocol. The high acidity of such adhesives compromises the chemical cure of the resin cement.
Metal-free, metal-ceramic or composite crown cementation
1.The full porcelain or metal-ceramic piece must be previously submitted to an aluminum oxide jet treatment by the laboratory in order to generate superficial microporosities and in the case of composite pieces, the preparation can also be done with fine burs. Prepare the ceramic piece (if its an etchable ceramic) applying fluoridric acid (e.g. Condac Porcelana, FGM) in its internal surface, according to the manufacturer´s indications. Rinse the piece abundantly with water to remove all the acid and dry it well.
2. Silanization of the piece: apply silane (Prosil, FGM), which is an uncolored solution, to the piece with the aid of a disposable microapplicator (Cavibrush, FGM) and wait for 1 minute. Dry the surface with slight jets of compressed air.
3. Select the appropriate color of dual-curing resin cement for cementation of the piece. Place the appropriate amount of Allcem on the spatulation plate and mix the pastes during 10 seconds with spatula or use the self-mixing tip which should be discarded after its use. Allcem is presented in a double-body syringe (base + catalyst pastes) and disposable self-mixing applicator tips, which make the mixture of the pastes in the adequate proportion without creating bubbles. Allcem is also presented in isolated syringes of base paste and catalyst paste.
4. Apply a fine layer of cement inside the piece or in the cavity covering the whole surface.
5. Carefully position the piece on the tooth allowing the cement to slowly flow from all the margins.
6. Keep the piece in position and start cleaning the excess cement approximately 2 minutes after beginning to mix the pastes. In case the excess cement is removed right after the positioning of the piece, the margins should be light cured to impede the effect of the layer inhibited by oxygen. In case of minimally retentive restorations, position the piece and keep it still during the cleaning process of the excess cement.
7. After removing the excess of cement, light cure each surface and margins for 40 seconds or allow chemical curing to occur for 10 minutes after mixing the pastes.
8. Proceed to finishing and polishing.
Note: after light curing of the restoration and its margins, remove the excess of cement with a scalpel. The use of rotatory instruments should be avoided before total curing of the cement which may occur in about 10 minutes after mixing the pastes, according to the local temperature.
9. Instruct the patient not to bite on the restored teeth at least for 15 minutes after cementation.
Cementation of intrarradicular posts
1. Prepare the endodontically treated tooth to receive the post.
2. Etch enamel and dentin for 15 seconds with phosphoric acid at 37% (Condac 37, FGM). Wash with water in abundance and make sure that all the acid has been removed. Dry the canal with the aid of absorbent paper cones and dry dentin/enamel with gentle air blasts leaving the tissue moist.
3. Apply the adhesive (e.g. Ambar, Ambar APS or Ambar Universal APS, FGM) with the aid of a disposable microaplicator (e.g. Cavibrush, FGM) with active friction inside the root canal for 15 seconds, covering all the walls with a thin layer of adhesive. Make sure there is no excess of adhesive inside the canal. Follow manufacturer´s instructions to evaporate the solvent with air and light cure it for the recommended time.
Note: ALLCEM is incompatible to self etching single-step adhesives (except Ambar Universal APS, FGM) because it can interfere on the chemical cure. Three-step adhesives (Acid+Primer+Adhesive) or two-step adhesives (Acid+Primer/Bond or Acid/Primer+Bond) are indicated.
4. Prepare the post (e.g. WhitePost, FGM) applying a thin layer of silane agent (e.g. Prosil, FGM) and wait for 1 minute for it to dry. Apply gentle air streams.
5. Apply the desired quantity of ALLCEM on a mixing pad and mix it with a spatula for about 10 seconds, or use the self mixing tip in case of the dual body syringe. Apply a thin layer of the mixed cement on the post´s surface and/or inside the prepared root canal with proper instruments.
6. Place the post inside the canal. Clean the exceeding cement in approximately 2-3 minutes after the mixture of the cement. Proceed to light curing for 40 seconds on the occlusal surface to keep the post stable in position. The use of translucent posts increases the amount of light that get inside the canal, what helps curing the cement.
7. Proceed to coronal reconstruction of the tooth.
Conservation and storage
Syringe must be immediately closed after each use.
Do not expose the product to high temperatures or intense light.
Store the product at temperatures between 5°C and 25°C / 41°F and 77°F.
Expiration date is printed on the packaging.
Warnings
All dental products should be used only by dentists, dental prosthetics and authorized laboratories. Do not use the product after it is expired. When disposing the product, follow the legislation of your country. Keep out of the reach of children.

This material was manufactured for dental use only and must be manipulated according to the instructions of use. The manufacturer is not responsible for any damage caused by misuse or improper manipulation of the product. The user is obligated to prove (before using the product and under his/her responsibility) that this product is compatible with the desired utilization, especially when such utilization is not indicated in the instructions of use. Data description does not constitute any type of warranty and therefore is not binding.

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