FGM

Condac37

PROFISSIONAL
PROFESSIONAL
PROFESSIONAL
PROFESSIONNEL
PROFESSIONELL
ПОФЕССИОНАЛЬНЫЙ
PROFESJONALNEGO

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DENTSCARE LTDA - Av. Edgar Nelson Meister, 474 - Distrito Industrial - 89219-501 - Joinville - SC - Brasil / Brazil. CNPJ / Tax ID: 05.106.945/0001-06 INDÚSTRIA BRASILEIRA / BRAZILIAN INDUSTRY - Responsável Técnico / Technical Authority: Friedrich Georg Mittelstädt CRQ: 13100147 - SC - ANVISA nº 80172310025 - Marca / Brand: FGM® / Fabricado no Brasil / Made in Brazil

EC REP NUNO FLORES - Al. Bonifacio Lazaro Lozano, 3 – Piso 0 – C 2780-125 Oeiras / Portugal · Tel: (351) 21 4439292 · ec.rep@fgm.ind.br

Serviço ao profissional {Brasil}
0800 644 6100
Technical support {Europe}
Servicio al profesional {Americas}
+55 47 3441 6100
www.fgm.ind.br | fgm@fgm.ind.br

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EN Instructions for use

37% Phophoric acid etchant for enamel and dentin

Professional use only

Read all the information in these instructions for use carefully before using the product. Keep them for consultation, at least until the product has been totally used or there is no more interaction of the product with its last patient.

Product description

Condac 37% is water-based gel containing Phosphoric acid at 37%, with low viscosity and thixotropic properties. It is blue, which makes it easier to be visualized and controlled during the application. It mixes well with water which allows it to be easily removed after the etching. The content of phosphoric acid (orthophosphoric) allows for each a good chemical treatment of the surface.

Commercial Availability

Check the commercially available packages in your country.

· 3 syringes with 2,5mL each and 3 applicator tips.

· 1 syringe with 12,5g and 5 applicator tips.

Basic composition

Phosphoric acid at 37%, thickener, pigment and deionized water.

Product indication

Etching of the enamel and dentin to improve the adhesion of the restorative materials to the tooth.

Precautions and contraindications

· Highly acid product.

· When handling the product, the professional should wear gloves and goggles.

· Check if the applicator tip is not blocked and that the syringe is functioning well before starting to apply the product.

· The patient should wear an apron so that there is no contact with the skin or clothes as well as goggles. The product should not be handled over or very near the patient.

· If there is contact with the skin and soft tissues of the oral cavity, wash the affected areas immediately with water.

· Avoid contact with the eyes. If it occurs by accident, wash them abundantly with water and, if necessary, seek medical assistance.

· Protect the surrounding dental surfaces with a polyester strip or the gingival barrier Top dam (FGM) in areas which are not to be conditioned.

· Not injectable.

· Do not reuse the syringe.

Side Effects

If the product remains on the dental surface for a longer time than that specified, it can cause excessive demineralization.

Instructions

1 - Prepare the area to be etched cleaning and drying it. Protect the cavity with appropriate material before the etching, if necessary.

2 - The gel must be applied for 15 seconds for both enamel and dentin.

3 - After application, the tooth must be rinsed abundantly with water in a way to completely remove the acid from the dental surface. The dentin must be dried although not dehydrated (remain humid), as recommendation of the adhesive system's protocol.

Conservation and storage

Maintain the product in its original packaging always well sealed. Store the product at 5ºC and 30ºC | 41ºF and 86°F. Do not freeze the product.

Expiration date

Printed on the product packaging.

Warnings

Do not use the product after its expiry date. In order to discard the product, follow the legislation of your country. Keep out of reach of children.

This material was only manufactured for dental use and must be handled according to its instructions for use. The manufacturer is not responsible for damage caused by other uses or incorrect handling. Besides this, the user is obliged to prove, before use and under his or her responsibility that this material is compatible with the use desired, mainly when this use is not indicated in the instructions for use. Data descriptions do not constitute any type of guarantee and so are not in way linked.

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