FGM CE

Desensibilize KF 2%

BR | ES | EN | DE

PROFISSIONAL
PROFESIONAL
PROFESSIONAL
PROFESSIONELL

DENTSCARE LTDA - Av. Edgar Nelson Meister, 474 - Distrito Industrial - 89219-501 - Joinville - SC - Brasil / Brazil, CNPJ / Tax ID: 05.106.945/0001-06 INDÚSTRIA BRASILEIRA / BRAZILIAN INDUSTRY - Responsável Técnico / Technical Authority: Friedrich Georg Mittelstädt CRQ: 13100147 - SC ANVISA nº 80172310018 - Marca / Brand: FGM® / Fabricado no Brasil / Made in Brazil

EC REP NUNO FLORES - Al. Bonifacio Lazaro Lozano, 3 – Piso 0 – C 2780-125 Oeiras / Portugal · Tel: (351) 21 4439292 · ec.rep@fgm.ind.br

Serviço ao profissional {Brasil}
0800 644 6100
Technical support {Europe}
Servicio al profesional {Américas}
+55 47 3441 6100
www.fgm.ind.br | fgm@fgm.ind.br

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Dessensibilizante

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EN Instructions for use

Desensitizer

Professional use only
Read all the information of these directions carefully before using the product. Keep them for consultation, at least until the product has been totally used or there is no more interaction of the product with your last patient.

Product description
Desensibilize KF 2% is a desensitizing gel of low viscosity with a base of potassium nitrate and sodium fluoride. It has neural action from the Potassium Nitrate and dentin tubules occlusive action from the Sodium Fluoride. It has a more powerful action on sensitivity and can interfere with the speed of the dental bleaching. It should only be used before and during the bleaching when there is great sensitivity and it is necessary to lower the penetration of the bleaching agent in the tooth in order to avoid the sensitivity. It is not advisable for the patient to use this product at home owing to the toxicity of the Sodium Fluoride and the ensuing risk if it is ingested.

Presentation form
Packaging containing a syringe with 2.5g of product and 1 application tip.

Basic composition
Active ingredients: Potassium Nitrate at 5% and Sodium Fluoride at 2%.

Inactive ingredients: Deionized Water, Glycerine, neutralizing agent and thickener.

Product indication
It is indicated in the treatment of teeth with moderate to high sensitivity to temperature variances and food (dentinal Hyperaesthesia) and in the treatment of sensitivity caused by dental bleaching.

Precautions and contraindications
· Avoid excessive ingestion of the product on purpose.
· Do not apply the product on gums in a bad state of health.
· Do not use it on patients with anaemia and/or any sensitivity to the product components.
· Discontinue using it if the patient has any undesirable reaction to the product components.
· Not injectable. Do not reuse the syringe.

Side Effects
Ingestion of the product in large quantities can cause nausea, intoxication by Fluorides and accentuate states of anaemia (effects of the Potassium Nitrate).

Instructions for use
1. Clean the surface of the teeth eliminating plaque, calculus and other deposits which can prevent the product from penetrating the teeth;
2. Apply the desensitizing gel uniformly on the teeth (use a brush or Cavibrush micro-applier) and wait for at least 10 minutes. Alternatively the gel can be applied with the help of the bleaching tray, executing the application as if it were the bleaching agent itself.
3. After the application time has ended, remove the gel from the teeth with cotton wool and plenty of water. Do not allow the patient to swallow the desensitizing agent.
4. Repeat the application of the desensitizing agent as required or, in the case of dental bleaching, following up the treatment time. The number of applications must not exceed once a day.

Note: After the application time of the desensitizing agent has ended and while it is still on the teeth, it is advisable to rub the product with the micro-applier, brush or felt disks at a low rotation in intermittent circular movements (low pressure). This procedure mechanically helps the desensitizing agent to penetrate the tooth and so enhances its performance.

Conservation and storage
Maintain the product in its original packaging always well sealed. Store the product at 5°C to 30°C | 41ºF to 86ºF. Do not freeze the product.

Expiration date
Printed on the product packaging.

Warnings
Do not use the product after expiry date. In order to discard the product, follow the legislation of your country. Keep out of reach of children.
This material was only manufactured for dental use and must be handled according to the instructions for use. The manufacturer is not responsible for damage caused by other uses or incorrect handling. Besides this, the user is obliged to prove, before use and under his or her responsibility, that this material is compatible with the use desired, mainly when this use is not indicated in the instructions for use. Data descriptions do not constitute any type of guarantee and so are not in way linked.

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