FGM

CE

Diamond R

BR | ES | EN | FR | DE

PROFISSIONAL | PROFESIONAL | PROFESSIONAL | PROFESSIONNEL | PROFESSIONELL

DENTSCARE LTDA - Av. Edgar Nelson Meister, 474 - Distrito Industrial - 89219-501 - Joinville - SC - Brasil / Brazil. CNPJ / Tax ID: 05.106.945/0001-06 INDUSTRIA BRASILEIRA / BRAZILIAN INDUSTRY - Responsável Técnico / Technical Authority: Friedrich Georg Mittelstädt CRQ: 13100147 - SC - ANVISA nº 80172310003 - Marca / Brand: FGM® / Fabricado no Brasil / Made in Brazil

EC REP NUNO FLORES - Al. Bonifacio Lazaro Lozano, 3 – Piso 0 – C 2780-125 Oeiras / Portugal · Tel: (351) 21 4439292 · ec.rep@fgm.ind.br

Serviço ao profissional {Brasil}
0800 644 6100
Technical support {Europe}
Servicio al profesional {Americas}
+55 47 3441 6100
www.fgm.ind.br | fgm@fgm.ind.br

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Pasta de polimento

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Polishing paste

Professional use only

Read all the information in this manual carefully before using the product. Keep the manual for later reference at least until the product is finished or until there is no more interaction between the product and its last patient.

Product description

The polishing paste DIAMOND R was developed in order to obtain polishing in resins. In its formulation it has carbowax, which acts as an excipient and lubricant for the polishing device avoiding generating too much heat and abrasives with very fine granulation. The paste is produced using atoxic ingredients, soluble in water, especially selected to help with the lubrication during the polishing (minimizing the generation of heat during the procedure) and facilitate its removal at the end of the treatment. The paste has a medium viscosity and is not affected by heat (it does not melt). It has an agreeable flavor of mint.

Presentation form

Packaging containing 1syringe of 4g.

Basic composition

Aluminum oxide, Carbowax, thickener, mint essence and water.

Product indication

It is indicated for use in the final polishing of resins in general

Precautions and contraindications

· Use low rotation and intermittent movements with a low application pressure of the felt on the tooth. This avoids excessive heating of the surface.

· Always use new felt for executing the polishing. Using the same felt disk with the abrasives of other pastes can jeopardize the success of the polishing.

Instructions for use

1.Prepare the area to be polished, leaving it smooth and in the end format desired. If necessary, use sanding disks (Diamond Pro Disks granulations G, M, F and/or XF) or medium and fine paste (Diamond ACI and ACII) for finishing prior to the polishing. This stage is of paramount importance, as obtaining a shine depends very much on the starting condition of the surface.

2.Deposit a small quantity of paste on a felt disk (Diamond or Diamond Flex Disks) or directly on the restoring material and proceed with the polishing.

3.Repeat procedure 2, if necessary, to obtain the effect desired.

Conservation and storage

Maintain the product in its original packaging always well sealed. Keep it out of direct sunlight. Store the product at 5°C to 30ºC / 41°F to 86ºF. Do not freeze the product.

Expiration date

Printed on the product packaging.

Warnings

Do not utilize the product after the expiry date. When disposing of the product, follow the legislation of your country. Keep out of children's reach.

This material was manufactured for dental use only and must be manipulated according to the manufacturer's instructions of use. The manufacturer is not responsible for any damage caused by misuse or improper manipulation of the product. The user is obligated to prove (before using the product and under his/her responsibility) that this product is compatible with the desired utilization, especially when such utilization is not indicated in the instructions of use. Data description does not constitute any type of warranty and therefore is not binding.

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