Case 2 | *Caso 2*

INITIAL/ *INICIAL*

FINAL/ *FINAL*

EN

36-year-old female patient

MAIN COMPLAINT

Esthetics of teeth 47 and 46 and bad breath.

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INITIAL CLINICAL / X-RAY EVALUATION

Amalgam restoration on tooth 47 and poor sealing on tooth 46 caused bad breath, because they allowed the accumulation of food debris with ease.

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TREATMENT PERFORMED

Removal of the restorations was performed and, after the full cleaning of the cavity, the protocol of absolute isolation of the teeth to be restored was carried out. It was decided to restore tooth 46, which had a larger depth and uneven bottom wall, with Opus Bulk Fill APS composite firstly in order to even the bottom wall and then Vittra APS (FGM) composite to stratify the restoration. Tooth 47 was restored entirely with Vittra APS (FGM) composite through increments, because it was a minor restoration.

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Clinical Tip:
Opus Bulk Fill Flow streamlines our day-to-day clinical practice because of its excellent fluidity and thixotropy.

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Fig. 8 - Opus Bulk Fill Flow APS (FGM) composite light-cured in the cavity.

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Fig. 9 - Placement of Vittra APS (FGM) composite in shade DA2 by the incremental technique.

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Clinical Tip:

To mask the possible translucency of the bulk fill composite used as a base, it is interesting to cover it with a higher opacity composite, as is the case with dentin composites.

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Fig. 11 - Tooth 46 completed.

Fig. 11 - Diente 46 finalizado.

Fig. 12 - Placement of dentin increment in tooth 47 (Vittra APS - FGM) in shade DA2.

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