EN Restorations of posterior teeth with bulk-fill composite combining esthetics with more agility.

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EN

44-year-old female patient

MAIN COMPLAINT
Pain in tooth 25 when eating sweet foods.

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INITIAL CLINICAL / X-RAY EVALUATION

Sensitivity to air jet upon clinical examination. From the interproximal x-ray image, a radio translucent image was observed indicating caries and infiltration surrounding the amalgam restoration.

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TREATMENT PERFORMED

Exchange of the amalgam restoration by restoration with bulk-fill composite (Opus Bulk Fill APS - FGM). The material was chosen due to practicality and esthetics. Besides not changing its shade after light curing, it allows us to work in bigger increments without losing its shape while being layered.

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Ambar Universal APS

Stability and longevity of color and adhesion of esthetic materials.

Self-etching light-curing adhesive system for enamel and dentin.

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Fig. 4 - After acid etching on enamel and application of the Ambar Universal APS (FGM) adhesive on enamel and dentin, the distal box was built with Opus Bulk Fill APS (FGM) composite in shade A1 and light-cured for 40 seconds.

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Fig. 5 - The matrix was then repositioned to the mesial area for building this box in a single increment.

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Fig. 6 - The proximal boxes were built to transform the cavity into a class I restoration.

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Fig. 7a and 7b - Single increment of Opus Bulk Fill APS (FGM) in shade A1, being "accommodated" in the cavity and applied up to the occlusal surface, thus restoring the anatomy of the tooth.

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Clinical Tip:

Opus Bulk Fill APS is perfect to be used in large increments, in addition to having low shrinkage stress, the APS system allows longer working time with the material.

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