EN Maximum safety and practicality allied with a bactericidal glass ionomer and low compression composite

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
20-year-old, female patient.

COMPLAINT

Patient reported sensitivity during chewing and ingestion of cold drinks in tooth 36.

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EN CLINICAL/RADIOGRAPHIC EXAM

An unsatisfactory composite restoration (class I of Black) was found in tooth 36 due to margin fracture, absence of functionally adequate anatomical shape, micro-infection in the tooth/restoration interface and color change. Both aesthetic and functional changes. The patient presented periodontal health with a wide range of keratinized tissue. At the probing examination, the depth of the sulcus was within the healthy standards. The radiographic exam suggested the presence of remnant carious tissue underlying the restoration.

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TREATMENT

Replacement of the unsatisfactory restoration was performed.


> *Opus Bulk Fill APS follows the worldwide trend for simplifying filling procedures, allowing restorations with shorter clinical time and smaller possibility of errors.*
>
> *Prof. Dr. Luis Fernando Morgan*

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STEP BY STEP | *PASO A PASO*

1a and 1b. Radiographic and clinical appearance of the case, showing unsatisfactory occlusal restoration.

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



2a to 2c. Removal of restorative material with diamond spherical tip number 1012 and conformation of cavity walls with carbide drill number 245, providing surrounding walls parallel to each other, rounded internal angles and clear cavosurface angle without a bevel.

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


3a and 3b. Note the unevenness of the pulp wall. In the deepest region, glass ionomer cement (Ion-Z) was applied as a liner, according to the principles of protection of the dentin pulp complex. The zinc present in the formulation of this product has a bactericidal action.
4a and 4b. Acid etching (Condac 37) only on enamel the margins for 30 seconds and subsequent removal with jets of water abundantly.

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

5a to 5c. With proper moisture control, application of the self-etching adhesive (Ambar Universal APS) to dentin and enamel and light curing for 20 seconds, leaving the cavity ready to receive restorative material.

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

6a and 6b. The cavity was filled with low-shrinkage stress composite (Opus Bulk Fill Flow APS) with a uniform distance of 2mm between the surface of composite and the cavosurface angle. Light curing was conducted for 40 seconds.

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7a to 7e. The cavity then received increments of composite (Opallis) with 2mm each in enamel shade. With this, the occlusal anatomy is re-established. The final increments were light cured for 20 seconds each and the final restoration for another 40 seconds.

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8a and 8b. Restoration completed after finishing and polishing (Diamond Flex and Diamond Excel). Note the interesting visual aspect of the shade effect obtained in the grooves and cusps by the combination of composites Opus Bulk Fill Flow APS and Opallis.

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LIGHT-CURING COMPOSITE
Opallis

RESINA FOTOCURABLE
Opallis

7TH GENERATION COLORLESS LIGHT-CURING ADHESIVE
Ambar Universal APS

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GLASS-IONOMER BASED RESTORATIVE MATERIAL
Ion-Z

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