FGM

CE 2797

OPUS BULK FILL FLOW APS

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EN Instructions for use

Flowable light-curing composite resin

Professional use only

Read carefully all the information on this instructions manual before using the product. Keep it for future reference at least up to the total consumption of the product or there is no further interaction with its last patient.

Product description

Opus Bulk Fill Flow APS is a bulk-fill flowable light-curing composite (fluid) recommended to be used as base and lining of restorations on posterior teeth. With approximately 68% of load filler, the composite shows enough mechanical strength to be used under direct restorations. The most superficial layer of the restoration (at least 2mm thick) should be restored with conventional methacrylate-based composite resins.

Forms of presentation

Refil: 1 syringe with NET WT. 0.07 OZ (2g) of resin and 10 tips.

Basic compostition

· Active ingredients: urethane dimethacrylate monomers, stabilizers, photoinitiating composition (APS) and coinitiator.
· Inactive ingredients: inorganic load filler of silanized silica dioxide, stabilizers and pigments.

Product indication:

· Base/lining underneath direct restorations.
· Small, non occlusal stress-bearing class I restorations according to minimally invasive filling therapy
· Pit and fissure sealant
· Tunnel-type preparation.
· Repair of enamel defects
· Bonding of tooth fragments.
· Repair in composite resin
· Non-carious cervical lesions.
· Planning of preparation walls.
· Core build-up

Precautions and contraindications

· For dentistry use only;
· Read carefully the recommendations and instructions for use;
· Use Opus Bulk Fill Flow APS only as informed in these instructions;
· Do not use Opus Bulk Fill Flow APS in patients with a history of allergy to composite resins that are acrylate based. In case of allergy, suspend the use and request the patient to look for medical assistance;
· In order to minimize the risk of irritation, regions of the cavity preparation too close to the pulp should be adequately protected (Eg.: application of calcium hydroxide);
· To avoid allergic reactions to methacrylate monomer, it is recommended to constantly wear gloves. If Opus Bulk Fill Flow APS gets in contact with the hands, wash them thoroughly with water;
· In case of eye contact, wash them thoroughly with water and contact a doctor immediately;
· When using a light-curing machine, eye protection should be worn.
· Opus Bulk Fill Flow APS is contraindicated when it is not possible to establish dry operative field or it is not possible to follow the described technique correctly;
·.Eugenol- and hydrogen peroxide containing materials should not be used in conjunction with this product because they may interfere with hardening and cause softening of the polymeric components of the material
· Opus Bulk Fill Flow APS should not be swallowed or inhaled. In case this occurs, search for medical help immediately;
·The enclosed intraoral tips are intended for a single use and must be disposed of after
only one use. Return the cap to the syringe after use.
. Contact with saliva and blood during composite placement may cause failure of the
restoration. Use of rubber dam or adequate isolation is recommended.

Side effects

· The product may cause some level of irritation if in contact with the dental pulp or too close to it. In cavities that are too deep, use a lining material with calcium hydroxide base;
· In cases of allergic reactions to the product, suspend its use. Opus Bulk Fill Flow APS is contraindicated for patients that are allergic or hypersensitive to methacrylate and related monomers and should not be used in those patients.
· Product may irritate the eyes and skin. Eye contact: irritation and possible corneal damage.
· Skin contact: irritation or possible allergic response. Reddish rashes may be seen on the skin.

Instructions for use

Note 1: Opus Bulk Fill Flow APS is recommended for use according to the adhesive technique associated to adhesives of one step, two steps or three steps.
Note 2: Opus Bulk Fill Flow APS should be used together with the conventional resins that are methacrylate-based, which should be applied as the last layer
1. Prepare the cavity and maintain the tooth away from contaminants, preferably with total dam isolation. If the cavity is too deep, use liners that are calcium-hydroxide based.
2. If the area involves proximal faces of the tooth (class II), use a metallic matrix and cervical wedge to allow adequate adaptation of the restoration.
3. The application of phosphoric acid and dental adhesive should follow the protocol of the respective products that are used. Opus Bulk Fill Flow APS is compatible with all the adhesive techniques. Once the hybrid layer is created, proceed to the next step
Application of the resin in the cavity
4.1 Opus Bulk Fill Flow APS should be applied to the teeth directly with the applicator tip that comes with the product. The tip should be threaded firmly to the syringe before usage, substituting the syringe lid;
4.2 After firmly coupling the tip, dispense a small amount of resin on a glass plate or mixing pad. That will guarantee that occasional air bubbles are removed from the syringe leaving the material ready for use;
4.3 Apply Opus Bulk Fill Flow APS in the cavity in increments of up to 4mm. Begin applying the resin on the deepest parts of the cavity, dispensing the material in a way that the applicator tip is on the same level or just above the layer that is being applied. Do not submerge the tip in the resin during application;
4.4 Light-cure each increment during 20 seconds using an adequate light-curing machine: minimum power of 450mW/cm2 with wave length of 400-500nm;
4.5 Fill the cavity until there is 2mm left or more for the occlusal surface. That last layer (the most superficial one) should be restored with a conventional composite resin.
5. Check occlusal contacts and proceed to finishing and polishing of the restoration.

Conservation and storage

Keep out of direct sunlight, exposure to light leads to premature curing.Store in a well ventilated place at temperatures between 5 °C and 27 °C. (41°F and 80,6°F). Allow material to reach room temperature prior to use. Protect from moisture. Do not freeze. Do

not use after expiration date.

Expiration date

Printed on the product packaging.

Warnings

Any dental product should only be used by dentists, prosthetic technicians and authorized dental laboratories. Do not use the product if it is expired. For disposal, follow the legislation of your country. Keep it out of the reach of children. Avoid contact of non-cured Opus Bulk Fill Flow APS with skin, mucosae and eyes. Before light-curing, the product may have a slight irritating effect and promote hypersensitivity due to methacrylates.

This material was manufactured for dental use only and must be handled according to the instructions for use. The manufacturer is not responsible for damages caused for other use or incorrect manipulation. Additionally, the user is obligated to verify, before the use and under his/her responsibility, if this material is compatible with the desirable usage, mainly in cases in which such usage is not indicated in this instructions for use. Data description do not constitute any type of warranty and, for this reason, are not binding.

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