FGM

CE 2797

OPUS BULK FILL CAPS

BR | ES | EN | FR | DE

PROFISSIONAL | PROFESIONAL | PROFESSIONAL | PROFESSIONNEL | PROFESSIONELL

DENTSCARE LTDA - Av. Edgar Nelson Meister, 474 - Distrito Industrial - 89219-501 - Joinville - SC - Brasil - CNPJ / Tax ID: 05.106.945/0001-06 INDÚSTRIA BRASILEIRA / BRAZILIAN INDUSTRY - Responsável Técnico / Technical Authority: Friedrich Georg Mittelstädt CRO: 13100147 - SC - ANVISA nº 80172310078 - Marca / Brand: FGM® / Fabricado no Brasil / Made in Brazil

EC REP NUNO FLORES - Al. Bonifacio Lazaro Lozano, 3 – Piso 0 – C 2780-125 Oeiras / Portugal · Tel: (351) 21 4439292 · ec.rep@fgm.ind.br

Serviço ao profissional (Brasil)
0800 644 6100
Technical support (Europe)
Servicio al profesional (Americas)
+55 47 3441 6100
Customer Service (USA)
1 888 827 7478 | cs@fgm.ind.br

www.fgm.ind.br | www.fgm-usa.com | fgm@fgm.ind.br

Irritante | Irritant | Reizend

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Resina composta fotopolimerizável

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EN Instructions for use

Light-curing composite resin

Professional use only
Read carefully all the information in this manual before using the product. Keep it for further reference, at least, until the product is thoroughly used or until there is no more interaction of the product with its last patient.

Product description
Opus Bulk Fill APS is a light-curing composite recommended for restorations performed in large increments. The composite has approximately 79% of load filler (in mass).
APS is the acronym for Advanced Polymerization System, and it consists of a combination of different photoinitiators that interact among each other.
The APS system is activated by any type of light-curing unit available on the market, even those with less power or narrow range of wavelength.

Presentation
Refill: 1 syringe with NET WT. 0.14 OZ (4g) of resin.
Caps: package containing 20 caps with NET WT. 0.007 OZ (0.2g) each.

Basic composition
· Active ingredients: uretane-dimetacrylic monomers, stabilizers, photoinitiating composition (APS) and co-initiators.
· Inactive ingredients: inorganic load of silanized silicon dioxide (silica), stabilizers and pigments.

Product indication
· Direct restorations in posterior teeth (permanent and deciduous) in increments of up to 5mm, including the occlusal surface;
·Class V restorations (cervical caries, root erosion, wedge-shaped defects)
· Reconstructive build-up
· Base and lining of direct restorations;
· Repair of small defects of the enamel;
· Repair of temporary acrylic and resin materials.

Precautions and contraindications
· This product is intended to be used only according to this instructions manual;
· Do not use the product if the patient is allergic to one of the substances present in the composition. In case of allergic reactions, discontinue the use of the product and, if needed, seek medical attention;
· Avoid the contact of uncured resin with the skin, mucosa and eyes because it may present an irritating effect due to the presence of methacrylates. It is recommended to wear all the Personal Protection Equipment (PPE) while manipulating the product;
· Avoid the use of eugenol containing materials, such as liners or temporary restoratives, as they hinder proper light curing of the product;
· When restoring deep cavities, lining material (such as calcium hydroxide) must be used on the pulp wall;
· The composite is sensitive to light, so make sure to keep the syringe properly closed when not in use. The resin content should be used right after being extruded from the syringe.;
· Use protective goggles when light curing the composite;
· The product must not be swallowed or aspirated. If so, seek medical attention;
· The product is contraindicated when a dry field cannot be stablished or proper technique cannot be followed;

Side effects
· The product may cause a certain degree of irritation when in contact with the dental pulp and/or when it gets too close to it. In cavities that are too deep, use a lining material that is calcium-hydroxide-based;
· In case of allergic reactions to the product, interrupt its use. Opus Bulk Fill APS is not recommended for patients that are allergic or hypersensitive to methacrylate and related monomer and should not be used in those patients.

Instructions for use
Note 1: Opus Bulk Fill APS is recommended for use according to the adhesive technique associated to single-step adhesives, two-step adhesives or three-step adhesives.
Note 2: Opus Bulk Fill APS may be used as a capping resin (occlusal surface) when associated to flowable bulk-fill resins.
1. Prepare the cavity and maintain the tooth free from contaminants preferably with rubber dam isolation. If the cavity is too deep, use linings that are calcium-hydroxide based.
2. If the cavity involves proximal faces of the tooth (class II), use metallic matrix and cervical wedge to allow for a good adaptation of the restoration;
3. The application of phosphoric acid and dental adhesive should follow the protocol of the respective products to be used. Opus Bulk Fill APS is compatible with all the adhesive techniques. Once the hybrid layer is created, proceed to the next step.
Application of resin to the cavity
4.1. Opus Bulk Fill APS should be applied to the tooth directly by using an appropriate spatula.
4.2 Apply Opus Bulk Fill APS in the cavity in increments up to 5mm deep. Begin applying the resin in the deepest areas of the cavity, condensing always against the dentin walls. It is possible to conclude the restoration with resin, including the occlusal surface.
4.3 Light cure during 40s when using a light-curing device of 450-1000 mW/cm² or during 30s when using a light-curing device of 1000-2000 mW/cm². It is suggested to distribute light-curing time pointing the light from different angles and to different surfaces (occlusal, buccal, lingual/palatal and proximal). The light-curing device must present wavelength of 400-500nm and its tip should be as close to the teeth/restoration as possible.
5. Verify occlusal contacts and proceed to the finishing and polishing of the restoration.

Single dose capsule:
1. Insert the capsule into the syringe delivery system properly, so that the capsule is ready to be used.
2. Remove the cover from the capsule and deliver the composite into the cavity by applying hand pressure to the gun/syringe.
3. After performing the restoration, discard the capsule.

Conservation and storage
The product should remain in its original packaging and be kept tightly closed;
Protect from direct sunlight;
Store the product at temperatures between 5°C-27°C / 41°F-80.6°F.
Do not freeze the product.

Expiration date
Printed on the product packaging.

Warnings
· Dental products should be used by authorized professionals only;
· Do not use the product after its expiration date;
· Follow your country's legislation when disposing the product;
· Keep out of the reach of children.

This material was manufactured for dental use only and must be manipulated according to the instructions for use. The manufacturer is not responsible for damages caused by other uses or incorrect manipulation. Additionally, the user is obligated to verify, before the use and under his/her responsibility, if this material is compatible with the desired utilization, mainly in cases in which such utilization is not indicated in this manual. Data description do not constitute any type of warranty and, for this reason, are not binding.

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