EN Installation of braces with self-adhesive cement.

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INITIAL CLINICAL / X-RAY EVALUATION

The patient had dental class II, mild dental rotation, unevenness and midline deviation. Tooth 26 was extruded, making it impossible to restore tooth 36 by reducing the occlusal space.

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TREATMENT PERFORMED

We chose orthodontic correction with a self-ligating metal active appliance. Buildups were installed to assist in the movement. The brackets were cemented with Orthocem (FGM) directly and without prior application of the adhesive. The resinous buildups were made with Orthobite (FGM) and, as they are colored, their removal is easy when necessary.

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Fig. 3 - Grabado ácido del esmalte.
Fig. 3 - Acid etching of enamel.

Orthocem

More comfort for the patient and greater profitability for your clinic.

Light-curing self-adhesive cement for bracket bonding.

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Fig. 7 - Brackets and arches installed.
Fig. 7 - Brackets y arcos instalados.

Ortho Bite

Resistance and comfort for your patient.

Light-curing resin for occlusal stops and band cementation.

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Source/fuente: Daniel Tocolini, Bruno Lippmann e Iduilton Grabowski.

Clinical Tip:

Top Comfort provides comfort to the patient as it protects soft tissues. In the picture, the product was applied to the hooks of the appliance to avoid injury to the mucosa.

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Light-curing orthodontic wax.
Cera ortodóntica fotocurable.

"The great clinical advantage of Orthocem (FGM) is the ease of installation of the brackets because of its excellent viscosity, which keeps the bracket in position before its light-curing.

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