FGM

CE

Top dam

PROFISSIONAL
PROFESIONAL
PROFESSIONAL
PROFESSIONNEL

BR | ES | EN | FR

DENTSCARE LTDA - Av. Edgar Nelson Meister, 474 - Distrito Industrial - 89219-501 - Joinville - SC - Brasil / Brazil. Autorização de funcionamento / Authorization of functioning: MS P5X44XY0XX28 - CNPJ / Tax ID: 05.106.945/0001-06 INDÚSTRIA BRASILEIRA / BRAZILIAN INDUSTRY - Responsável Técnico / Technical Authority: Friedrich Georg Mittelstädt CRQ: 13100147 - SC Registro na ANVISA / Registration at ANVISA nº 80172310023 - Marca / Brand: FGM® / Fabricado no Brasil / Made in Brazil

EC REP NUNO FLORES - Al. Bonifacio Lazaro Lozano, 3 – Piso 0 – C 2780-125 Oeiras / Portugal · Tel: (351) 21 4439292 · ec.rep@fgm.ind.br

Serviço ao profissional {Brasil}
0800 644 6100
Technical support {Europe}
Servicio al profesional {Americas}
+55 47 3441 6100
www.fgm.ind.br | fgm@fgm.ind.br

Irritante / Irritant

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Barreira gengival fotopolimerizável

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Light cured resin barrier

Professional use only

Read all the information in this manual carefully before using the product. Keep the manual for later reference at least until the product is finished or until there is no more interaction between the product and its last patient.

Product description

Top dam is a light cured resin which replaces the rubber dam (complete isolation) in some cases, where the isolation related to the gingiva is sufficient. The product is sufficiently viscous and thixotropic to allow easy application and control over the place to be protected. It has a bluish colour in order to make a

contrast with thegingiva and the teeth and so facilitate its localization. After covering the area which you wish to protect, just polymerize the product with a blue light and it becomes rigid, insoluble and slightly adhering to the gingiva. At the end of the treatment the polymerized product can be easily detached and removed.

Presentation form

Packaging containing a syringe of 2g and three disposable tips for applying the product.

Basic composition

HEMA, monomers of urethane di-methacrylate, inert charge, pigments and photoinitiators.

Product indication

It is indicated as a replacement of the rubber sheet to protect the gingival tissue in cases where there is a risk of irritation by contact with products used in the clinical procedure, such as for example dental bleaching with Hydrogen Peroxide and to help with the complete isolation for better sealing of a cavity preparation or endodontical treatment.

Precautions and contraindications

· Avoid applying the product on gingival tissue in poor conditions, people who have undergone recent periodontal surgery or which contains a previously installed irritation.

· The incomplete polymerization of the product or excessive time can lead to a slight gingival irritation in some cases. If the patient has a stronger allergic reaction to the product, it suffices to interrupt using the product.

· Not injectable. Do not reuse the syringe.

Side Effects

The product contains methacrylate monomers which can cause allergic reactions or slight irritations in patients who are sensitive to these substances.

During the process of polymerizing the product a slight sensation of heat or burning may be felt on the gingiva owing to the reaction of the polymerization and contraction of the product.

Instructions

1.Dry the gingiva and the teeth in the area which you wish to protect.

2.Apply a fine layer of the product (thicicness of 1 - 1,5mm) on the gingiva and the interdental papilla in a way that they are completely covered and there is a good field of protection for the procedure to be applied.

3. Light cured 20 to 30 seconds. It is advisable to proceed with the polymerization soon after applying the resin in order to prevent any possible irritation.

4.Verify that the product covers all the gingiva surface and that it is completely polymerized. In order to see the rear part of the papillae and gingiva, use a clinical mirror looking from incisal to cervical. Avoid testing it by forcing it or trying to lift it as this can cause infiltration of the irritating substance.

5.In order to remove the product just move it forwards. It should come loose easily as a whole or in large pieces.

Conservation and storage

Keep the product in its original packaging always well sealed. Keep it out of direct sunlight. Store at 5°C to 35ºC / 41°F to 95ºF. Do not freeze the product.

Expiration date

Printed on the product packaging.

Warnings

Do not utilize the product after the expiry date. When disposing of the product, follow the legislation of your country. Keep out of children's reach.

This material was manufactured for dental use only and must be manipulated according to the manufacturer's instructions of use. The manufacturer is not responsible for any damage caused by misuse or improper manipulation of the product. The user is obligated to prove (before using the product and under his/her responsibility) that this product is compatible with the desired utilization, especially when such utilization is not indicated in the instructions of use. Data description does not constitute any type of warranty and therefore is not binding.

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