FINAL

**The result after whitening shows the success of the mixed technique.**

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INITIAL | INICIAL

## CASE HIGHLIGHTS | DESTAQUES DEL CASO

*Microabrasion can be performed manually or mechanically with the aid of a rubber cup at low rotation speed. Each application should last 10-15 seconds (per tooth) and then the teeth should be washed copiously.*

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*As microabrasion slightly increases the roughness of the enamel, it is important to complete the session with polishing.*

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*In this case, 5 applications were performed in the same session. Note the more homogeneous color of the teeth.*

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EN

# **Microabrasion** combined with **tooth whitening**

AUTHORS: Claudio Sato and Adriano Sapata

*Female patient, 29.*

## COMPLAINT

Yellow and stained teeth.

## INITIAL APPEARANCE

Hypoplastic stains very evident on the upper incisors, characteristic of dental fluorosis.

## CLINICAL/X-RAY EXAM

White stains on enamel were mainly observed in the areas of enamel dominance: the incisal edges of the maxillary anterior teeth. Apparent on the surface, with well-defined edges, present bilaterally, suggesting dental fluorosis.

## TREATMENT

Microabrasion of enamel was associated with tooth whitening. When microabrasion is used, the removal (even micrometrically) of a superficial layer of enamel makes the color of the dentin even more apparent. This fact suggests that microabrasion should be combined with tooth whitening for superior results. The active tips used in the treatment may be rotary (rubber cups or cones) or manual (spatulas, wooden wedges). The choice of the size of these active tips always depends on the size of the stain, as tips much larger than the stain may provide unnecessary wear of healthy and stainless areas.

> " Microabrasion is an effective, practical and painless procedure. In cases of light staining, in association with take-home whitening it has a great effect. "
>
> *DR. CLAUDIO SATO AND DR. ADRIANO SAPATA*

## FGM MATERIALS USED

*STAIN REMOVER BY MICROABRASION*
***Whiteness RM***

*DIAMOND POLISHING PASTE*
***Diamond Excel***

*TAKE-HOME DENTAL BLEACHING GEL (CARBAMIDE PEROXIDE)*
***Whiteness Perfect 10%***

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