Result achieved: drastic reduction in staining, shade 1M1 for canines and 0.5M1 for the other teeth.

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INITIAL | INICIAL

CASE HIGHLIGHTS | DESTAQUES DEL CASO

Before microabrasion, the whitening was carried out for just one week for mitigation of stains and to increase the value of the teeth.

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The fluorotic enamel has a lower hardness than healthy enamel, and thus has lower resistance to wear during abrasion.

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The application of Desensibilize KF2% desensitizing agent made the process more comfortable, since the product contains potassium nitrate, which interferes with the depolarization of nerve fibers, thereby inhibiting pain stimuli and reducing sensitivity. Furthermore, it also contains sodium fluoride, which acts as a remineralizing agent. Desensibilize KF2% desensitizing agent was also used directly by the patient before the whitening with trays.

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EN

Minimally invasive techniques for dental fluorosis: **microabrasion and whitening**

AUTHORS: Alexander Nishida and Carlos Francci

Female patient, 20.

COMPLAINT

Teeth with brown stains and many white stains.

INITIAL APPEARANCE

Initial appearance of fluorosis stains. The brownish stains on the anterior teeth were causing embarrassment to the patient and were her main complaint.

CLINICAL/X-RAY EXAM

Stains on teeth were classified as grade 5 fluorosis in accordance with the TF index (Thylstrup & Fejerskov).

TREATMENT

After prophylaxis/scaling and shade taking, supervised take-home whitening, it was carried out with 10% carbamide peroxide for 4 hours daily (only one syringe) aimed at smoothing the stains and whitening the substrate. By decreasing the difference between the color of the stains and the color of the teeth, less microabrasion would be needed. Microabrasion was divided into two sessions: the first consisting of 10 applications on each upper arch and the second session covering the lower arch. Desensitizing agent was used for greater comfort to the patient. When the microabrasion step was completed, take-home whitening was performed again with the same concentration.

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The use of Diamond Pro sandpaper disks improved the finishing and polishing of proximal and cervical areas due to the availability of two diameters, 8 and 12mm. Treatment should always begin with the more abrasive disk, which has a more saturated blue color, moving on to the other less abrasive disks, with the blue color becoming lighter.

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" Microabrasion and dental whitening allow resolution of the aesthetic problems in a minimally invasive way, with a high degree of satisfaction for the patients and for the professionals, and at a low cost. "

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DR. ALEXANDER NISHIDA AND DR. CARLOS FRANCCI

FGM MATERIALS USED

STAIN REMOVER BY MICROABRASION
Whiteness RM

DESENSITIZING AGENTS
Desensibilize KF2% and KF0.2%

DIAMOND POLISHING PASTE
Diamond excel

PLATES FOR TRAYS
Whiteness Plates for Custom Dental Trays 1mm

TAKE-HOME DENTAL BLEACHING GEL (CARBAMIDE PEROXIDE)
Whiteness Perfect 10%

SODIUM FLUORIDE AT 1.23% FOR TOPICAL USE
Flúor Care

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