FINAL

After the microabrasion session, the surface stains were removed.

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INITIAL | INICIAL

CASE HIGHLIGHTS | DESTAQUES DEL CASO

Abrasion can be done mechanically, with rubber cups, which allows for control of the pressure on the teeth with greater precision. In this case, 13 applications were performed, lasting 5 seconds each, in one session.

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After microabrasion, fluoride is applied topically to aid in remineralization. This fluoride must be applied for 1 minute and afterward, the patient must refrain from eating for 30 minutes.

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EN

Aesthetic resolution **after orthodontic treatment** with the technique of **microabrasion** of the enamel

AUTHORS: Nathália Nunes Vidal, Juliana Yassue Barbosa da Silva and Sheila de Carvalho Stroppa

Male patient, 15.

COMPLAINT

Stains on the teeth after removal of braces.

INITIAL APPEARANCE

White stains apparent, especially on the upper anterior teeth.

CLINICAL/X-RAY EXAM

Teeth 13, 12, 11, 21, 22, 23, and 24 showed active demineralization on vestibular enamel located around the site where the orthodontic brackets were placed.

TREATMENT

After guidance on oral hygiene and remineralization of the active white stains on enamel with 4 sessions of topical application of fluoride sealant and after rebalancing of the oral environment, the microabrasion technique was carried out with Whiteness RM (FGM) stain remover on the respective teeth for aesthetic restoration of the patient.

> "Microabrasion removes the stains painlessly and definitively, in a practical and quick procedure."
>
> *DR. SHEILA DE CARVALHO STROPPA*

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FGM MATERIALS USED

STAIN REMOVER BY MICROABRASION
Whiteness RM

SODIUM FLUORIDE AT 1.23% FOR TOPICAL USE
Flúor Care

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